

Protocols for CRISPR sgRNA

For genome engineering in cultured cells and animal embryo injection

1. Reconstitution of lyophilized sgRNA

- a) Dissolve lyophilized guide RNA in nuclease-free water
- b) Recommended stock solution: 2-5 ug/ul
- c) Store sgRNA stock solutions at - 80°C

2. Validation of sgRNA activity

- a) Set up reaction mix:

Components	Amount
Cas9 protein	500ng
sgRNA	250ng
Targeting substrate	80~200ng
10 x NEB Buffer 3.1	1µl
10 x BSA	1µl
DW	to 10µl

- b) Incubate the reaction mixture at 37°C for 1 hr
- c) Add 1 ul of RNase and incubate for 15 min at 37°C
- d) Add 1 ul of STOP solution (30% glycerol, 1.2% SDS, 250 mM EDTA (pH 8.0)) to the reaction mixture and incubate for 15 min at 37°C
- e) Analyze on 2% agarose gel

3. Gene Knockout Cell Establishment

The amount of the reagents given in the protocol below are for one well of a 24-well plate. For other reaction formats, scale the amounts of reagents up or down accordingly.

- a) Add Cas9 RNP complex (0.5 µg of Cas9 Nuclease and 250 ng of CRISPR sgRNA) to 50 µl Opti-MEM I Reduced Serum Medium
- b) In a separate tube, dilute the transfection reagent by adding 4 µl of Lipofectamine 2000 transfection reagent to 50 µl of Opti-MEM I Reduced Serum Medium
- c) Mix gently and incubate for 5 min at RT

- d) Add the diluted transfection reagent to the tube containing Cas9 protein/sgRNA RNP complexes and mix gently
- e) Incubate at room temperature for 20 minutes to allow the formation of Cas9/sgRNA complex
- f) Add the Cas9/sgRNA complex to cells to be transfected (e.g. 1×10^5 NIH3T3)
- g) Swirl the plates gently to allow the mixing of the transfection mixture with the medium
- h) Incubate the plate at 37°C in a humidified CO₂ incubator in a cell culture incubator for 2~3 days
- i) Assay samples to determine the genome editing efficiency by T7E1 assay or targeted deep sequencing

4) General guidelines for the establishment of a KO animal by direct injection of CRISPR RNP complex into a one-cell mouse embryo

- a) Make a mixture of Cas9 protein and sgRNA:
Cas9 at 20 ng/ul and sgRNA 10 ng/ul (1 ul each) are the recommended concentrations for mouse embryos.
- b) Incubate for 15 min at 37°C
- c) Inject diluted solution into one-cell mouse embryos